



# COLD HANDLING CABINET



Spacious cold handling cabinet made of stainless steel for semen processing at working temperature +5°C. Timer for automatic defrosting cycle. Cover made of acrylic glass. Installed on heavy duty rolls.

Available in 2 sizes:

	working area (cm - W x D)	outside dimensions (cm - W x D x H)	230 V	115 V
large	200 x 90	216 x 114 x 127	REF. : 14335/0403	REF. : 14335/0400
standard	200 x 67	216 x 91 x 127	REF. : 14335/0401	REF. : 14335/0398

*Shorter, custom made versions are available.*

## Technical features

- Optionally with hydraulic height adjustment
- Height of refrigerated area: 500 mm (340 mm below acrylic cover in the front)
- Digital temperature controller
- 4 electric connectors in the working area
- 50 mm openings on each side for cables
- Power supply: 230 V/50 Hz or 115 V/60 Hz
- Power consumption: approx. 450 W

## Your benefits

- Homogeneous temperature distribution throughout the whole cabinet
- Precise temperature control
- Extra space for loading freezing racks
- Comfortable working height
- Robust and maintenance free

